

Part List - 2026

P/N	Description	Notes to Consider
	Tube Mount	
CM-11-0021-2B	CM Tube Mount Assembly	<i>Includes the Mounting Collar</i>
CP-11-SR21-4A	SourceRay Tube Mount Assembly	<i>Includes the Mounting Collar</i>
CP-11-MX21-3A	MinXRay Tube Mount Assembly	<i>Includes the Mounting Collar</i>
CM-12-0018-1D	Mounting Collar	<i>Sold as assembly only</i>
CM-12-MX18-1A	MinXRay Mounting Collar	<i>Sold as assembly only</i>
CP-12-0115-1A	Aluminum Filter, 0.25mm	<i>Package of 5</i>
CM-12-0123-1B	Mounting Spacer	<i>Package of 4</i>
CM-12-0162-1A	Tube Mount Retaining Plate	<i>Sold as assembly only</i>
	Housing	
CM-11-0128-1B	Housing Assembly, Beige	<i>Includes: Lead, Output Window, Bezel, Tape Measure, Overlay</i>
CP-11-0126-1A	Housing Assembly, White	<i>Includes: Lead, Output Window, Bezel, Tape Measure, Overlay</i>
CM-11-0153-1A	Housing Assembly, White White Overlay w/ Grey Lettering	<i>Includes: Lead, Output Window, Bezel, Tape Measure, White Overlay</i>
CP-11-SR05-1F	Housing Assembly, SourceRay	<i>Includes: Lead, Output Window, Bezel, Tape Measure, Overlay</i>
CP-11-MX05-1H	Housing Assembly, MinXRay (w/ Skin Guard)	<i>Includes: Lead, Output Window, Skin Guard, Angle Indicator, Tape Measure, All Overlays</i>

CM-11-0004-1B	CM Tape Measure	<i>Includes Mounting Adhesive</i>
CP-11-MX04-1B	MinXRay Tape Measure	<i>Includes Mounting Adhesive</i>
CP-11-SR04-1B	SourceRay Tape Measure	<i>Includes Mounting Adhesive</i>
CP-11-0080-1C	Output Window Assembly	<i>Includes Mounting Adhesive</i>
CM-11-0120-1D	Bezel Assembly	<i>No longer Available</i>
CP-11-0007-1D	Skin Guard Assembly	<i>Includes Mounting Adhesive</i>
	Top Cover	
CP-12-0006-1B	Top Cover, Beige	
CM-11-0006-1B	Top Cover Assembly, Beige	<i>Includes the Switch and Wire Leads</i>
CP-12-0172-1B	Top Cover, White	
CP-11-0127-1B	Top Cover Assembly, White	<i>Includes the Switch and Wire Leads</i>
CP-12-MX06-1B	Top Cover, Black, MinXRay	
CP-11-MX06-1B	Top Cover Assembly, Black, MinXRay	<i>Includes the Switch and Wire Leads</i>
CM-11-0145-1A	Lamp/Laser Switch Assembly	<i>Includes Wire Leads</i>
	Electrical Cord	
CP-13-0045-1A	34" Collimator Power Cord No Connector, 3 Wire	
CM-11-0139-1A	Adapter, 3 Pin Connector, Standard	<i>If included on a new collimator PO, we will replace the 9 pin connector at no charge</i>

CM-11-0143-1A	Adapter, 3 Pin Connector, Locking	<i>If included on a new collimator PO, we will replace the 9 pin connector at no charge</i>
CM-11-0046-1B	34" Collimator Power Cord 3 Pin Connector, Locking	
CM-11-0134-1A	34" Collimator Power Cord 9 Pin Connector	
CM-11-0047-1B	25' Interface Power Cable 3 Pin Connector, Locking	
CM-11-0121-1A	25' Interface Power Cable, 9 Pin Connector	
CM-11-0135-1A	50' Interface Power Cable 9 Pin Connector	
CM-11-0136-1A	65' Interface Power Cable 9 Pin Connector	
CML-13-0001	24VDC Power Supply	<i>Includes 110V AC Wall Outlet Power Cord with Spade Terminals and Wiring Instructions</i>
	Turn Knob	
CM-11-0138-1A	Turn Knob Assembly	<i>Does Not Include Indicators Package of 2</i>
CM-11-0142-1A	Turn Knob Assembly, Tapered Fit	<i>Does Not Include Indicators, Package of 2 For collimators with S/N less than A04586</i>
CM-11-0119-1A	Indicator Assembly	<i>Fits All Collimare Knobs Package of 2</i>
CP-12-0036-1B	Indicator (Machined)	<i>No Longer Available Replace with CM-11-0119-1A</i>
CP-17-0036-1A	Indicator (As Molded)	<i>No Longer Available Replace with CM-11-0119-1A</i>
CM-11-0002-1A	Turn Knob Assembly, Aluminum Core	<i>No Longer Available Replace with CM-11-0138-1A</i>
CM-11-0060-1A	Turn Knob Assembly Black Cap, Gray Plastic Core	<i>No Longer Available Replace with CM-11-0138-1A</i>
CP-99-0034-1A	Gray Plastic Turn Knob Package of 2	<i>No Longer Available Replace with CM-11-0142-1A</i>

	Overlay	
CM-16-0096-1A	LED Face Plate Overlay	<i>No Longer Available Replace with CM-16-0101-1A</i>
CM-16-0097-1A	LED Face Plate Overlay Metric	
CM-16-0100-1B	LED Face Plate Overlay Metric, White	
CM-16-0082-1B	LED Face Plate Overlay 60in SID	
CM-16-0101-1A	Smart' LED Face Plate Overlay	<i>Gray Background w/ White Lettering Standard Overlay for CMLs</i>
CM-16-0102-1A	'Smart' LED Face Plate Overlay White	<i>White Background w/ Grey Lettering</i>
CP-16-MX82-1A	MinXRay Overlay Scale Side	
CP-16-MX83-1A	MinXRay Overlay Face Plate	
CP-16-MX84-1A	MinXRay Overlay Blank Side	
CP-16-SR90-1B	LED Face Plate Overlay Source-Ray	<i>No Longer Available Replace with CP-16-SR91-1A</i>
CP-16-SR91-1A	'Smart' LED Face Plate Overlay Source-Ray	
	Accessory	
CM-11-0123-1A	Internal Heavy Weight Kit	<i>Heavy Weight Collimator Only P/N: CML-150-HV01-C-B2</i>
CM-10-1001-1B	Accessory Rail Kit	<i>Includes Hardware</i>
CM-10-1002-1B	DAP Meter Rail Kit	<i>Includes Hardware</i>
CM-12-0166-1A	Siemens Tube Mount Adapter	<i>Includes Hardware</i>

	Chassis	
Repair	Internal Shutter/Cable Repair	<i>Collimator has to be returned to factory</i>
Loaner	Loaner Collimator	<i>Loaner Collimator For Use During Collimator Factory Repair</i>
CP-11-0011-1E	Mounting Lug Assembly	<i>Includes: Adjustment Screw, Mounting Lug, Mounting Slide</i>
CP-12-0167-1A	1/4-20 Double Acting Adjustment Screw	<i>Part of Mounting Lug Assembly Package of 4</i>
CP-11-0122-1A	Drive Cam	<i>Sold as Drive Shaft Assembly Only</i>
CP-11-0046-1A	Transverse Drive Shaft	<i>Includes Drive Cam</i>
CP-11-0047-1A	Longitudinal Drive Shaft	<i>Includes Drive Cam</i>
	Mirror	
CP-12-0042-1A	Light Field Mirror	
CP-14-0009-1A	Mirror Retaining Spring Clip	<i>Package of 2</i>
	Laser	
CP-12-0047-2A	Laser Module Housing	<i>No Longer Available Replace with CP-13-0010-1A</i>
CP-13-0010-1A	Laser Module Assembly Internal Adjustment	<i>For collimators with S/N less than A08323 Includes Laser Module Housing</i>
CP-11-0124-1A	Front Laser Assembly External Adjustment	<i>For collimators with S/N greater than A08322 Front to Back Laser Line</i>
CP-11-0125-1A	Side Laser Assembly External Adjustment	<i>For collimators with S/N greater than A08322 Left to Right Laser Line</i>
	Circuit Board	

CM-11-0673-1E	Standard LED Board Assembly	<i>No Longer Available Replace with CM-11-A674-1A</i>
CM-11-0674-1A	High Heat LED Board Assembly	<i>For collimators with S/Ns greater than A14406 Replaces an existing High Heat LED Board</i>
CM-11-A674-1A	High Heat LED Board Field Upgrade	<i>Replaces Standard LED Board (S/N less than A14407) Includes Replacement Instructions</i>
CP-11-SR673-1A	Standard LED Board Assembly, Source Ray	<i>No Longer Available Replace with CP-11-SR674-1A</i>
CP-11-SR674-1A	High Heat LED Board, Source Ray Field Upgrade	<i>Source Ray Models, Includes Mounting Brackets Replaces Standard LED Board (S/N less than A14407)</i>
CP-13-0001-1A	Halogen Circuit Board	<i>CM/CP Halogen Models</i>
CP-13-0077-1A	LED Power Conditioning Board	<i>CML 150/125/VET Models For use with a Standard LED Board</i>
CP-13-0066-1A	Connector. 6 Screw, 10 amp	<i>For LED Board</i>
CP-13-0080-1A	Connector. 8 Screw, 10 amp	<i>For LED Board</i>
	Halogen Units	
CP-13-0030-1B	100W 12V lamp	<i>MinXRay Models</i>
CP-13-0031-1A	100W 24V lamp	<i>SourceRay Models</i>
CP-13-0032-1A	150W 24V lamp	<i>CM Models</i>
CP-11-MX21-2B	MinXRay Tube Mount Assembly, Floating Mount	<i>Includes the Mounting Collar Smooth Ring, No Posts</i>
CP-11-SR21-3A	SourceRay Tube Mount Assembly, Floating Mount	<i>Includes the Mounting Collar Smooth Ring, No Posts</i>
CP-12-0018-2E	Mounting Collar, Floating Mount	<i>Sold as Assembly Only Smooth Ring, No Posts</i>
CP-12-MX18-3B	MinXRay Mounting Collar, Floating Mount	<i>Sold as Assembly Only Smooth Ring, No Posts</i>

CP-13-0020-1A	Thermal Switch	
CM-12-0122-1B	Heat Sink Cover	
CM-11-0032-1B	Halogen Aperture Plate Assembly	<i>Includes Lamp Base</i>
CP-13-0040-1A	Lamp Base	
	Collimators	
CML-150-0001-C	Manual, 150kVp, High Heat 'Smart' LED, 24VDC, 24VAC	<i>Crosshair Lasers, Swivel Tube Mount, High Heat 'Smart' LED Control with programmable on time up to 12 minutes, 25' Interface Cable with remote switch capability, Installation and Instruction Kit</i>
CML-150-HV01-C	Manual, Heavy Weight, 150kVp, High Heat 'Smart' LED, 24VDC, 24VAC	<i>Crosshair Lasers, Swivel Tube Mount, High Heat 'Smart' LED Control with programmable on time up to 12 minutes, 25' Interface Cable with remote switch capability, Installation and Instruction Kit. Direct Replacement for Heavier Collimators (Progeny, Eureka, Ralco, etc) - Approx 18lbs</i>
CML-125-0001-C	Manual, 125kVp, High Heat 'Smart' LED, 24VDC, 24VAC	<i>Crosshair Lasers, Swivel Tube Mount, 'Smart' LED Control with programmable on time up to 12 minutes, 25' Interface Cable with remote switch capability, Installation and Instruction Kit</i>
CML-VET-0001-NC	Manual, Non-Certified "VET", High Heat 'Smart' LED, 24VDC, 24VAC	<i>Crosshair Lasers, Swivel Tube Mount, 'Smart' LED Control with programmable on time up to 12 minutes, 25' Interface Cable with remote switch capability, Installation and Instruction Kit</i>
Collimare Halogen Exchange	-\$350 off of the price of a LED Collimator	<i>The \$350 discount will be applied to the PO. If we do not receive back the halogen unit within 30 days from the ship date of the new collimator then \$350 will be billed.</i>